

DnaK ATPase-BD E.Coli

Description: Recombinant DnaK Substrate Binding Domain produced in E.Coli is a single, non-glycosylated polypeptide chain containing 384 amino acids and having a molecular mass of 48.1 kDa.

Catalog #: HYP5-017

Synonyms: HSP-70, HSP70, DnaK, Chaperone protein dnaK, Heat shock protein 70, Heat shock 70 kDa protein, groP, grpF, seg, b0014, JW0013.

For research use only.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MGKIIGIDLG TTNSCVAIMD GTTPRVLENA EGDRTTPSII
AYTQDGETLV GQPAKRQAVTNPQNTLFAIK RLIGRRFQDE EVQRDVSIMP FKIIAADNGD
AWVEVKGQKM APPQISAEVLKMKKTAEDY LGEPVTEAVI TVPAYFNDAQ RQATKDAGRI
AGLEVKRIIN EPTAAALAYGLDKGTGNRTI AVYDLGGGTF DISIIEIDEV DGEKTFEVL
TNGDTHLGGE DFDSR

Purity: Greater than 95.0% as determined by (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.

Formulation:

The DnaK protein contains 25mM Tris-HCl, pH7.5, 100mM NaCl, 5mM DTT and 10% Glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

DnaK, originally identified for its DNA replication by bacteriophage I in E. coli is the bacterial HSP-70 chaperone. This protein is involved in the folding and assembly of newly synthesized polypeptide chains and in preventing the aggregation of stress-denatured proteins. DnaK (amino acids 1-384) is N-terminal ATPase domain and ATP bound to the ATPase domain induces a conformational change in the substrate binding domain (residues 385-638). The protein coding region of the ATPase domain of DnaK (amino acids 1-384) was amplified by PCR and cloned into an E. coli expression vector. The ATPase domain of DnaK was purified to apparent homogeneity by using conventional column chromatography techniques.

To place an order, please [Click HERE](#).