

CAPN1

Description: CAPN1 consists of an 80-kDa large subunit and a 30 kDa small subunit. CAPN1 was purified by sequential chromatography through DEAE-Sepharose, A1.5m Bio-Gel, and Phenyl-Sepharose CL-4B columns.

Catalog #: PRPS-370

For research use only.

Synonyms: Calpain-1 catalytic subunit, EC 3.4.22.52, Calpain-1 large subunit, Calcium-activated neutral proteinase 1, Calpain mu-type, muCAPN, Micromolar-calpain, Cell proliferation-inducing gene 30 protein, CANP 1, CAPN1, CANPL1, PIG30, CANP, muCL, CANP1.

Source: Human Erythrocytes.

Physical Appearance: Sterile Filtered colorless solution.

Purity: Greater than 90% as determined by SDS-PAGE.

Formulation:

50mM imidazole-HCl, 100mM NaCl, 5mM EGTA, 1mM DTT and 10% sucrose.

Stability:

CAPN1 although stable at 10°C for 1 week, should be stored below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

This protein can be used for immunoblots, absorption experiments in immunohistochemistry, radioimmunoassay and intracellular injection. For adsorption we suggest the following procedure:

A- Dilute 1

Introduction:

Calpains activity is attributed to two main isoforms: μ -calpain and m-calpain, which are ubiquitously expressed proteases implicated in cellular migration, cell cycle progression, degenerative processes and cell death. These heterodimeric enzymes are composed of distinct catalytic subunits, encoded by Capn1 (μ -calpain) or Capn2 (m-calpain), and a common regulatory subunit encoded by Capn4. CAPN1 is a calcium-regulated non-lysosomal thiol-protease which catalyzes limited proteolysis of the substrates involved in cytoskeletal remodeling and signal transduction. CAPN1 is activated by micromolar concentrations of calcium and inhibited by calpastatin.

To place an order, please [Click HERE](#).