

ENO2

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 47kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human ENO2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NSE; HEL-S-279;

Catalog #: A0022

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2026

Isotype: IgG

Swiss Prot: P09104

Purity: Affinity purification

For research use only.

Background:

This gene encodes one of the three enolase isoenzymes found in mammals. This isoenzyme, a homodimer, is found in mature neurons and cells of neuronal origin. A switch from alpha enolase to gamma enolase occurs in neural tissue during development in rats and primates.

To place an order, please [Click HERE](#).