

RPS6KA1

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC IP FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 ICC 1:50 - 1:200 IP 1:20 - 1:50
FC 1:20 - 1:50

Calculated MW: 83kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human RPS6KA1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fpq

Synonym:

HU-1;MAPKAPK1A;RSK;RSK1;Ribosomal protein S6 kinase alpha-1 antibody; Ribosomal S6 kinase 1 antibody;

Catalog #: A0135

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6195

Isotype: IgG

Swiss Prot: Q15418

Purity: Affinity purification

For research use only.

Background:

RSK 1 (p90-RSK 1) is a member of a family of 90 kDa ribosomal protein S6 kinases, which includes RSK 1, RSK 2 and RSK 3. These are broadly expressed serine/threonine protein kinases activated in response to mitogenic stimuli, including extracellular signal-regulated protein kinases Erk1 and Erk2. p90-RSK 1 is activated by MAPK in vitro and in vivo via phosphorylation. Active RSKs appear to play a major role in transcriptional regulation by translocating to the nucleus and phosphorylating c-Fos and CREB.

To place an order, please [Click HERE](#).