

ARAF

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 68kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant Protein of human ARAF

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PKS2; A-RAF; ARAF1; RAFA1;

Catalog #: A0180

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 369

Isotype: IgG

Swiss Prot: P10398

Purity: Affinity purification

For research use only.

Background:

This proto-oncogene belongs to the RAF subfamily of the Ser/Thr protein kinase family, and maybe involved in cell growth and development. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).