

TNFRSF10D

Reactivity: Human

Tested applications: WB FC

Recommended Dilution: WB 1:500 - 1:2000 FC 1:20 - 1:50

Calculated MW: 42kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human TNFRSF10D

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bp

Synonym:

DCR2; CD264; TRUNDD; TRAILR4; TRAIL-R4;

Catalog #: A0184

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8793

Isotype: IgG

Swiss Prot: Q9UBN6

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor contains an extracellular TRAIL-binding domain, a transmembrane domain, and a truncated cytoplasmic death domain. This receptor does not induce apoptosis, and has been shown to play an inhibitory role in TRAIL-induced cell apoptosis.

To place an order, please [Click HERE](#).