

ATG5

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 32kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ATG5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bch

Synonym:

ATG5;APG5;APG5-LIKE;APG5L;ASP;hAPG5 ;

Catalog #: A0203

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9474

Isotype: IgG

Swiss Prot: Q9H1Y0

Purity: Affinity purification

For research use only.

Background:

Autophagy is a catabolic process for the autophagosomic-lysosomal degradation of bulk cytoplasmic contents (1,2). Autophagy is generally activated by conditions of nutrient deprivation but has also been associated with a number of physiological processes including development, differentiation, neurodegeneration, infection, and cancer (3). The molecular machinery of autophagy was largely discovered in yeast and referred to as autophagy-related (Atg) genes. Formation of the autophagosome involves a ubiquitin-like conjugation system in which Atg12 is covalently bound to Atg5 and targeted to autophagosome vesicles (4-6). This conjugation reaction is mediated by the ubiquitin E1-like enzyme Atg7 and the E2-like enzyme Atg10 (7,8).

To place an order, please [Click HERE](#).