

BID

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 22kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human BID

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bi

Synonym:

BID;FP497;MGC15319;MGC42355 ;

Catalog #: A0210

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 637

Isotype: IgG

Swiss Prot: P55957

Purity: Affinity purification

For research use only.

Background:

The BH3 domain-only protein, BID, a death agonist member of the Bcl-2/Bcl-xL family (1), is localized in the cytosolic fraction of cells as an inactive precursor (2,3). Its active form is generated upon proteolytic cleavage by caspase-8 in the Fas signaling pathway. Cleaved BID translocates to mitochondria and induces cytochrome c release and mitochondrial damage (2-5). Thus, BID relays an apoptotic signal from the cell surface to mitochondria. However, the precise molecular mechanism for the translocation of the cleaved BID, and for the subsequent release of cytochrome c during apoptosis, is still unclear.

To place an order, please [Click HERE](#).