

BVES

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 41kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant Protein of human BVES

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BVES ; Popeye protein 1; POPDC1; POP1; HBVES

Catalog #: A0213

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11149

Isotype: IgG

Swiss Prot: Q8NE79

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the POP family of proteins containing three putative transmembrane domains. This gene is expressed in cardiac and skeletal muscle and may play an important role in development of these tissues. The mouse ortholog may be involved in the regeneration of adult skeletal muscle and may act as a cell adhesion molecule in coronary vasculogenesis. Three transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Dec 2010]

To place an order, please [Click HERE](#).