

CASP10

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 59kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CASP10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

CASP10; MCH4; Caspase 10; FLICE2; CASP 10; ICE like apoptotic protease 4

Catalog #: A0216

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 843

Isotype: IgG

Swiss Prot: Q92851

Purity: Affinity purification

For research use only.

Background:

Caspase-10 is a DED (death effector domain)-containing caspase and functions as an initiator caspase in Fas/TNF induced apoptosis (1). Four isoforms of caspase-10 have been identified: caspase-10a (Mch4), caspase-10b (FLICE2), caspase-10c and caspase-10d. They have the same prodomain but different mature large and small subdomains (2-4). Upon death ligand-receptor binding, caspase-10 is coupled to the multimeric Fas/TNF receptor complex via DED/FADD adaptor interaction (1-4). This complex processes procaspase-10 into a large active fragment and a small fragment. Cleaved caspase-10 further processes other caspase members, including caspase-3 and caspase-7, to initiate a caspase cascade, leading to apoptosis (3-6).

To place an order, please [Click HERE](#).