

DDIT3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human DDIT3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Catalog #: A0221

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1649

Isotype: IgG

Purity: Affinity purification

For research use only.

To place an order, please [Click HERE](#).