

BMP2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 44kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human BMP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

prt

Synonym:

bmp2 ; BMP 2A; BMP 2; BMP2A; bone morphogenetic protein 2; Bone morphogenetic protein 2A

Catalog #: A0231

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 650

Isotype: IgG

Swiss Prot: P12643

Purity: Affinity purification

For research use only.

Background:

Bone morphogenetic proteins (BMPs) were first identified as molecules that can induce ectopic bone and cartilage formation (1,2). BMPs are synthesized as precursor proteins that are processed by cleavage to produce mature proteins. BMPs initiate signaling by binding to a receptor complex containing type I and type II serine/threonine receptor kinases that then phosphorylate Smad (mainly Smad1, 5 and 8), resulting in the translocation of Smad to the nucleus. BMP was also reported to activate MAPK pathways in some systems (3,4).

To place an order, please [Click HERE](#).