

VCAM1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:3000

Calculated MW: 81kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human VCAM1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VCAM1 ; L1CAM; VCAM 1; V CAM 1; DKFZp779G2333; CD106

Catalog #: A0279

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7412

Isotype: IgG

Swiss Prot: P19320

Purity: Affinity purification

For research use only.

Background:

Vascular cell adhesion protein 1 (VCAM1) is important in cell-cell recognition. It has been reported to function in leukocyte-endothelial cell adhesion, interact with the beta-1 integrin VLA4 on leukocytes, and mediate both adhesion and signal transduction. There are two isoforms; Isoform 1, a single-pass type I membrane protein, and isoform 2, a lipid-anchor (GPI-anchor). VCAM1 is expressed on inflamed vascular endothelium, as well as on macrophage-like and dendritic cell types in both normal and inflamed tissue. It is also expressed in the bone marrow (1).

To place an order, please [Click HERE](#).