

TREM1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:26kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human TREM1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TREM1;TREM-1 ;

Catalog #:A0292

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:54210

Isotype:IgG

Swiss Prot:Q9NP99

Purity:Affinity purification

For research use only.

Background:

This gene encodes a receptor belonging to the Ig superfamily that is expressed on myeloid cells. This protein amplifies neutrophil and monocyte-mediated inflammatory responses triggered by bacterial and fungal infections by stimulating release of pro-inflammatory chemokines and cytokines, as well as increased surface expression of cell activation markers. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.[provided by RefSeq, Jun 2011]

To place an order, please [Click HERE](#).