

XPO1

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 123kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human XPO1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

exp1; CRM1; exp1;

Catalog #: A0299

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7514

Isotype: IgG

Swiss Prot: O14980

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene mediates leucine-rich nuclear export signal (NES)-dependent protein transport. Exportin 1 specifically inhibits the nuclear export of Rev and U snRNAs. It is involved in the control of several cellular processes by controlling the localization of cyclin B, MPAK, and MAPKAP kinase 2. This protein also regulates NFAT and AP-1.

To place an order, please [Click HERE](#).