

MMP2

Reactivity: Human Mouse Rat

Tested applications: IHC

Recommended Dilution: IHC 1:50 - 1:100

Calculated MW: 72kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MMP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MMP2; CLG4; CLG4A; MMP-II; MONA; TBE-1 ;

Catalog #: A0317

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4313

Isotype: IgG

Swiss Prot: P08253

Purity: Affinity purification

For research use only.

Background:

The matrix metalloproteinases (MMPs) are a family of proteases that target many extracellular proteins including other proteases, growth factors, cell surface receptors and adhesion molecules (1). Among the family members, MMP-2, MMP-3, MMP-7 and MMP-9 have been characterized as important factors for normal tissue remodeling during embryonic development, wound healing, tumor invasion, angiogenesis, carcinogenesis and apoptosis (2-4). MMP activity correlates with cancer development (2). One mechanism of MMP regulation is transcriptional (5). Once synthesized, MMP exists as a latent proenzyme. Maximum MMP activity requires proteolytic cleavage to generate active MMPs by releasing the inhibitory propeptide domain from the full length protein (5).

To place an order, please [Click HERE](#).