

CD79A

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 25kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CD79A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD79A; IGA; MB1;

Catalog #: A0331

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 973

Isotype: IgG

Swiss Prot: P11912

Purity: Affinity purification

For research use only.

Background:

The B lymphocyte antigen receptor is a multimeric complex that includes the antigen-specific component, surface immunoglobulin (Ig). Surface Ig non-covalently associates with two other proteins, Ig-alpha and Ig-beta, which are necessary for expression and function of the B-cell antigen receptor. This gene encodes the Ig-alpha protein of the B-cell antigen component. Alternatively spliced transcript variants encoding different isoforms have been described.

To place an order, please [Click HERE](#).