

CCNB1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 48kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human CCNB1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

d

Synonym:

CCNB;

Catalog #: A0381

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 891

Isotype: IgG

Swiss Prot: P14635

Purity: Affinity purification

For research use only.

Background:

Cyclin B1 is a regulatory protein involved in mitosis. The gene product complexes with p34(cdc2) to form the maturation-promoting factor (MPF). Two alternative transcripts have been found, a constitutively expressed transcript and a cell cycle-regulated transcript, that is expressed predominantly during G2/M phase of the cell cycle. The different transcripts result from the use of alternate transcription initiation sites. The antibody is specific to CCNB1.

To place an order, please [Click HERE](#).