

KRT18

Reactivity: Human Mouse

Tested applications: WB IHC IF FC

Recommended Dilution: WB 1:1000 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200 FC 1:20 - 1:50

Calculated MW: 48kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human KRT18

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bj

Synonym:

PIG46;CYK18;K18; CK-18;

Catalog #: A0388

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3875

Isotype: IgG

Swiss Prot: P05783

Purity: Affinity purification

For research use only.

Background:

Keratins are a large family of proteins that form the intermediate filament cytoskeleton of epithelial cells. Keratin expression is highly regulated, tissue specific, and varies according to cell-state. Type I keratins consist of acidic, low molecular weight proteins with MW ranging from 40 kDa (KRT19) to 64 kDa (KRT9). Type 2 keratins consist of basic or neutral, high molecular weight proteins with MW from 52 kDa (KRT8) to 67 kDa (KRT18). KRT18, also named as CYK18, PIG46 and K18, is the most commonly found members of the intermediate filament gene family. KRT18 is a type I keratin which is found primarily in non squamous epithelia and is present in a majority of adenocarcinomas and ductal carcinomas but not in squamous cell carcinomas.

To place an order, please [Click HERE](#).