

KRT18

Reactivity: Human Rat

Tested applications: WB ICC

Recommended Dilution: WB 1:1000 - 1:4000 ICC 1:50 - 1:200

Calculated MW: 45kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human KRT18

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bj

Synonym:

K18; CYK18;

Catalog #: A0389

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 3875

Isotype: IgG

Swiss Prot: P05783

Purity: Affinity purification

For research use only.

Background:

KRT18 encodes the type I intermediate filament chain keratin 18. Keratin 18, together with its filament partner keratin 8, are perhaps the most commonly found members of the intermediate filament gene family. They are expressed in single layer epithelial tissues of the body. Mutations in this gene have been linked to cryptogenic cirrhosis. Two transcript variants encoding the same protein have been found for this gene.

To place an order, please [Click HERE](#).