

EIF4E

Reactivity: Human Mouse

Tested applications: WB FC

Calculated MW: 217

Observed MW: Refer to figures

Immunogen:

Synthetic Peptide

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A0467

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 1977

Isotype: IgG

Swiss Prot: P06730

Purity: Affinity purification

For research use only.