

DFFA

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 37kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human DFFA

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

DFF1; ICAD; DFF-45;

Catalog #: A0507

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1676

Isotype: IgG

Swiss Prot: O00273

Purity: Affinity purification

For research use only.

Background:

Apoptosis is accompanied by shrinkage and fragmentation of the cells and nuclei and degradation of the chromosomal DNA into nucleosomal units. DNA fragmentation factor (DFF), heterodimer of 40-kDa (DFFB) and 45-kDa (DFFA) subunits, is one of the major endonucleases responsible for internucleosomal DNA cleavage during apoptosis. DFFA is the substrate for caspase-3 and triggers DNA fragmentation during apoptosis. DFF activated after cleavage of DFFA by caspase-3 followed by its dissociation from DFFB. DFFB has been found to trigger both DNA fragmentation and chromatin condensation during apoptosis. DFFA appears to be a 45-kDa protein despite the predicated MW of 37 kDa.

To place an order, please [Click HERE](#).