

IRF2

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:39kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human IRF2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

IRF-2;

Catalog #:A0508

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3660

Isotype:IgG

Swiss Prot:P14316

Purity:Affinity purification

For research use only.

Background:

IRF2 encodes interferon regulatory factor 2, a member of the interferon regulatory transcription factor (IRF) family. IRF2 competitively inhibits the IRF1-mediated transcriptional activation of interferons alpha and beta, and presumably other genes that employ IRF1 for transcription activation. However, IRF2 also functions as a transcriptional activator of histone H4.

To place an order, please [Click HERE](#).