

RALA

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 24kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human RALA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MGC48949; RAL

Catalog #: A0541

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5898

Isotype: IgG

Swiss Prot: P11233

Purity: Affinity purification

For research use only.

Background:

The product of this gene belongs to the small GTPase superfamily, Ras family of proteins.

GTP-binding proteins mediate the transmembrane signaling initiated by the occupancy of certain cell surface receptors. This gene encodes a low molecular mass ras-like GTP-binding protein that shares about 50% similarity with other ras proteins. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).