

NEDD4

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 149kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NEDD4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

hlt

Synonym:

RPF1; NEDD4-1;

To place an order, please [Click HERE](#).

Catalog #: A0552

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4734

Isotype: IgG

Swiss Prot: P46934

Purity: Affinity purification

For research use only.