

SKP1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 19kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human SKP1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

d

Synonym:

EMC19; MGC34403; OCP-II; OCP2; SKP1A; TCEB1L; p19A

Catalog #: A0554

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6500

Isotype: IgG

Swiss Prot: P63208

Purity: Affinity purification

For research use only.

Background:

SKP1(S-phase kinase-associated protein 1) is also named as EMC19, OCP2, SKP1A, TCEB1L, belongs to the SKP1 family. It is an essential component of the SCF (SKP1-CUL1-F-box protein) ubiquitin ligase complex, which mediates the ubiquitination of proteins involved in cell cycle progression, signal transduction and transcription. This antibody is specific to SKP1.

To place an order, please [Click HERE](#).