

MAPK12

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW: 42kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human MAPK12

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

io

Synonym:

ERK3; ERK6; P38GAMMA; PRKM12; SAPK-3; SAPK3

Catalog #: A0559

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6300

Isotype: IgG

Swiss Prot: P53778

Purity: Affinity purification

For research use only.

Background:

Activation of members of the mitogen-activated protein kinase family is a major mechanism for transduction of extracellular signals. Stress-activated protein kinases are one subclass of MAP kinases. The protein encoded by this gene functions as a signal transducer during differentiation of myoblasts to myotubes. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).