

SOX2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 34kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SOX2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

mg

Synonym:

SOX2; ANOP3; MCOPS3; MGC2413 ;

Catalog #: A0561

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6657

Isotype: IgG

Swiss Prot: P48431

Purity: Affinity purification

For research use only.

Background:

Embryonic stem cells are derived from the inner cell mass of the blastocyst and are unique in their pluripotent capacity and potential for self-renewal. Sox2 is one of a set of transcription factors that are crucial for the maintenance of pluripotency (1). Sox2, Oct-4, and Nanog cooperate in this network (1-3), and siRNA knockdown of either Sox2 or Oct-4 results in loss of pluripotency (4,5). Chromatin immunoprecipitation experiments have shown that Sox2 and Oct-4 bind to thousands of gene regulatory sites, highlighting the importance of these transcription factors in early embryonic development (6,7). It has recently been shown that Sox2 is amplified in lung and esophageal squamous cell tumors (8).

To place an order, please [Click HERE](#).