

CENPA

Reactivity: Human Rat

Tested applications: WB IHC ICC IP FC

Calculated MW: 140

Observed MW: Refer to figures

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A0570

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1058

Isotype: IgG

Swiss Prot: P49450

Purity: Affinity purification

For research use only.