

SNAI2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:20 - 1:200

Calculated MW: 30kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human SNAI2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

SLUG; WS2D; SLUGH1; SNAIL2;

Catalog #: A0572

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6591

Isotype: IgG

Swiss Prot: O43623

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the Snail family of C2H2-type zinc finger transcription factors.

The encoded protein acts as a transcriptional repressor that binds to E-box motifs and is also likely to repress E-cadherin transcription in breast carcinoma. This protein is involved in epithelial-mesenchymal transitions and has antiapoptotic activity. Mutations in this gene may be associated with sporadic cases of neural tube defects.

To place an order, please [Click HERE](#).