

AXIN1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW: 96kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human AXIN1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AXIN1; AXIN; Axin-1; Axis inhibition protein 1; PPP1R49;

Catalog #: A0583

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8312

Isotype: IgG

Swiss Prot: O15169

Purity: Affinity purification

For research use only.

Background:

This gene encodes a cytoplasmic protein which contains a regulation of G-protein signaling (RGS) domain and a dishevelled and axin (DIX) domain. The encoded protein interacts with adenomatosis polyposis coli, catenin beta-1, glycogen synthase kinase 3 beta, protein phosphate 2, and itself. This protein functions as a negative regulator of the wingless-type MMTV integration site family, member 1 (WNT) signaling pathway and can induce apoptosis. The crystal structure of a portion of this protein, alone and in a complex with other proteins, has been resolved. Mutations in this gene have been associated with hepatocellular carcinoma, hepatoblastomas, ovarian endometrioid adenocarcinomas, and medullablastomas. Two transcript variants encoding distinct isoforms have been identified for this gene.

To place an order, please [Click HERE](#).