

CAMK2A

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:20 - 1:200

Calculated MW: 54kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human CAMK2A

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CAMKA; KIAA0968

Catalog #: A0585

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 815

Isotype: IgG

Swiss Prot: Q9UQM7

Purity: Affinity purification

For research use only.

Background:

The product of this gene belongs to the serine/threonine protein kinases family, and to the Ca(2+)/calmodulin-dependent protein kinases subfamily. Calcium signaling is crucial for several aspects of plasticity at glutamatergic synapses. This calcium calmodulin-dependent protein kinase is composed of four different chains: alpha, beta, gamma, and delta. The alpha chain encoded by this gene is required for hippocampal long-term potentiation (LTP) and spatial learning. In addition to its calcium-calmodulin (CaM)-dependent activity, this protein can undergo autophosphorylation, resulting in CaM-independent activity. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Nov 2008]

To place an order, please [Click HERE](#).