

YES1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW: 61kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human YES1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Yes; c-yes; HsT441; P61-YES

Catalog #: A0628

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7525

Isotype: IgG

Swiss Prot: P07947

Purity: Affinity purification

For research use only.

Background:

This gene is the cellular homolog of the Yamaguchi sarcoma virus oncogene. The encoded protein has tyrosine kinase activity and belongs to the src family of proteins. This gene lies in close proximity to thymidylate synthase gene on chromosome 18, and a corresponding pseudogene has been found on chromosome 22.

To place an order, please [Click HERE](#).