

SSB

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:100 - 1:200

Calculated MW: 47kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SSB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

La; LARP3

Catalog #: A0630

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6741

Isotype: IgG

Swiss Prot: P05455

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is involved in diverse aspects of RNA metabolism, including binding and protecting poly(U) termini of nascent RNA polymerase III transcripts from exonuclease digestion, processing 5' and 3' ends of pre-tRNA precursors, acting as an RNA chaperone, and binding viral RNAs associated with hepatitis C virus. Autoantibodies reacting with this protein are found in the sera of patients with Sjogren syndrome and systemic lupus erythematosus.

To place an order, please [Click HERE](#).