

CDC73

Reactivity: Human Mouse

Tested applications: WB IHC IF FC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:200 IF 1:20 - 1:50 FC 1:20 - 1:50

Calculated MW: 61kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human CDC73

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C1orf28; FLJ23316; HPT-JT; HRPT2

Catalog #: A0636

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79577

Isotype: IgG

Swiss Prot: Q6P1J9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a tumor suppressor that is involved in transcriptional and post-transcriptional control pathways. The protein is a component of the the PAF protein complex, which associates with the RNA polymerase II subunit POLR2A and with a histone methyltransferase complex. This protein appears to facilitate the association of 3' mRNA processing factors with actively-transcribed chromatin. Mutations in this gene have been linked to hyperparathyroidism-jaw tumor syndrome, familial isolated hyperparathyroidism, and parathyroid carcinoma. [provided by RefSeq, Jul 2009]

To place an order, please [Click HERE](#).