

GSK3A

Reactivity: Human Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 51kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GSK3A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

qt

Synonym:

DKFZp686D0638

Catalog #: A0645

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2931

Isotype: IgG

Swiss Prot: P49840

Purity: Affinity purification

For research use only.

Background:

This gene encodes a multifunctional Ser/Thr protein kinase that is implicated in the control of several regulatory proteins including glycogen synthase, and transcription factors, such as JUN. It also plays a role in the WNT and PI3K signaling pathways, as well as regulates the production of beta-amyloid peptides associated with Alzheimer's disease. [provided by RefSeq, Oct 2011]

To place an order, please [Click HERE](#).