

TOP2A

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:50

Calculated MW: 174kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human TOP2A

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bd

Synonym:

TOP2; TP2A

Catalog #: A0726

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7153

Isotype: IgG

Swiss Prot: P11388

Purity: Affinity purification

For research use only.

Background:

TOP2A, also named as TOP2, belongs to the type II topoisomerase family. Topoisomerase II makes double-strand breaks. TOP2 control the topological states of DNA by transient breakage and subsequent rejoining of DNA strands. The antibody is specific to the gene TOP2A.

To place an order, please [Click HERE](#).