

ERBB4

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 147kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human ERBB4

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ERBB4;HER4;MGC138404;p180erbB4 ;

Catalog #: A0750

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2066

Isotype: IgG

Swiss Prot: Q15303

Purity: Affinity purification

For research use only.

Background:

The HER/ErbB receptor tyrosine kinase family has been implicated in normal development, cardiac function and cancer (1,2). HER4/ErbB4, like other family members, has four ectodomains, a single transmembrane domain and a cytoplasmic tail containing the active tyrosine kinase domain (3). By binding to neuregulins and/or EGF family ligands, ErbB4 forms either a homodimer or heterodimer with other ErbB family members, which results in receptor activation and signaling (3). ErbB4 is ubiquitously expressed with the highest expression occurring in brain and heart. The expression of ErbB4 in breast cancer, pediatric brain cancer and other types of carcinomas has been reported suggesting that ErbB4 expression is involved in both normal tissue development and carcinogenesis (3).

To place an order, please [Click HERE](#).