

CD34

Reactivity: Human

Tested applications: WB IHC

Calculated MW: 385

Observed MW: Refer to figures

Immunogen:

Synthetic Peptide

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A0762

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 947

Isotype: IgG

Swiss Prot: P28906

Purity: Affinity purification

For research use only.