

VAV1

Reactivity: Human

Tested applications: WB ICC IP

Recommended Dilution: WB 1:500 - 1:2000 ICC 1:20 - 1:50 IP 1:20 - 1:50

Calculated MW: 95kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human VAV1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VAV;

Catalog #: A0771

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7409

Isotype: IgG

Swiss Prot: P15498

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the VAV gene family. The VAV proteins are guanine nucleotide exchange factors (GEFs) for Rho family GTPases that activate pathways leading to actin cytoskeletal rearrangements and transcriptional alterations. The encoded protein is important in hematopoiesis, playing a role in T-cell and B-cell development and activation. The encoded protein has been identified as the specific binding partner of Nef proteins from HIV-1. Coexpression and binding of these partners initiates profound morphological changes, cytoskeletal rearrangements and the JNK/SAPK signaling cascade, leading to increased levels of viral transcription and replication. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

To place an order, please [Click HERE](#).