

EEF1A1

Reactivity: Human

Tested applications: WB IHC IF FC

Calculated MW: 462

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human EEF1A1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A0831

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1915

Isotype: IgG

Swiss Prot: P68104

Purity: Affinity purification

For research use only.