

BIRC3

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:67kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human BIRC3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1mg/ml

Synonym:

AIP1; API2; CIAP2; HAIP1; HIAP1; MALT2; MIHC; RNF49

Catalog #:A0833

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:330

Isotype:IgG

Swiss Prot:Q13489

Purity:Affinity purification

For research use only.

Background:

Inhibitor of Apoptosis (cIAP) proteins are key intrinsic regulators of apoptosis induced by a variety of triggers. cIAPs overexpression has been shown to suppress apoptosis induced by a variety of stimuli. cIAP2 encodes a 60 kDa protein, originally identified as a molecule recruited to the TNFR complex, and expressed in many different cells and tissues.

To place an order, please [Click HERE](#).