

PRKAB1

Reactivity: Human Mouse Rat

Tested applications: WB IHC ICC IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:20 - 1:50 IP 1:50 - 1:200

Calculated MW: 38kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PRKAB1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

c

Synonym:

AMPK; HAMPKb;

Catalog #: A0846

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 5564

Isotype: IgG

Swiss Prot: Q9Y478

Purity: Affinity purification

For research use only.

Background:

Moesin (for membrane-organizing extension spike protein) is a member of the ERM family which includes ezrin and radixin. ERM proteins appear to function as cross-linkers between plasma membranes and actin-based cytoskeletons. Moesin is localized to filopodia and other membranous protrusions that are important for cell-cell recognition and signaling and for cell movement.

To place an order, please [Click HERE](#).