

FOXO1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 80kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human FOXO1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bk

Synonym:

FKH1;FKHR;FOXO1A;

Catalog #: A0885

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 2308

Isotype: IgG

Swiss Prot: Q12778

Purity: Affinity purification

For research use only.

Background:

This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in myogenic growth and differentiation. Translocation of this gene with PAX3 has been associated with alveolar rhabdomyosarcoma.

To place an order, please [Click HERE](#).