

AICDA

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 24kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human AICDA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

AID; ARP2; CDA2; HIGM2

Catalog #: A0899

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57379

Isotype: IgG

Swiss Prot: Q9GZX7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a RNA-editing deaminase that is a member of the cytidine deaminase family. The protein is involved in somatic hypermutation, gene conversion, and class-switch recombination of immunoglobulin genes. Defects in this gene are the cause of autosomal recessive hyper-IgM immunodeficiency syndrome type 2 (HIGM2). [provided by RefSeq, Feb 2009]

To place an order, please [Click HERE](#).