

PARK2

Reactivity: Mouse Rat

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:200

Calculated MW: 52kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PARK2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1mg/ml

Synonym:

PARK2;AR-JP;LPRS2;PDJ;PRKN;Parkin

Catalog #: A0968

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5071

Isotype: IgG

Swiss Prot: O60260

Purity: Affinity purification

For research use only.

Background:

Parkin is a protein of 465 amino acids with an amino-terminal ubiquitin domain and a carboxy-terminal RING-box (1). In the case of autosomal recessive juvenile Parkinsonism (AR-JP), deletions have been found in the gene on chromosome 6 encoding the protein Parkin (2).

To place an order, please [Click HERE](#).