

SIRT7

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 45kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SIRT7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

SIRT7;MGC126840;MGC126842;SIR2L7 ;

Catalog #: A0979

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51547

Isotype: IgG

Swiss Prot: Q9NRC8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the sirtuin family of proteins, homologs to the yeast Sir2 protein. Members of the sirtuin family are characterized by a sirtuin core domain and grouped into four classes. The functions of human sirtuins have not yet been determined; however, yeast sirtuin proteins are known to regulate epigenetic gene silencing and suppress recombination of rDNA. Studies suggest that the human sirtuins may function as intracellular regulatory proteins with mono-ADP-ribosyltransferase activity. The protein encoded by this gene is included in class IV of the sirtuin family. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).