

NCK1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 43kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NCK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NCK1;MGC12668;NCK;NCKalpha ;

Catalog #: A0999

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4690

Isotype: IgG

Swiss Prot: P16333

Purity: Affinity purification

For research use only.

Background:

NCK1 (also known as NCK or NCK) is a broadly expressed oncogenic adapter protein consisting of three SH3 domains and one SH2 domain (1-3). NCK1 becomes phosphorylated upon activation of variety of cell surface receptors and is involved in actin cytoskeletal organization induced by many stimuli (4-6). NCK2 (also known as NCK), a homolog of NCK1, has an overlapping expression pattern and redundant functions with NCK1 (7).

To place an order, please [Click HERE](#).