

PGP

Reactivity: Human

Tested applications: WB IHC FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 FC 1:20 - 1:50

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human PGP

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A10002

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 283871

Isotype: IgG

Swiss Prot: A6NDG6

Purity: Affinity purification

For research use only.