

EV71 3C

Reactivity:Enterovirus

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:5000

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human EV71 3C

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A10003

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Isotype:IgG

Purity:Affinity purification

For research use only.

To place an order, please [Click HERE](#).